

OCTs & Fundus Photography

Explore what they are & how they allow us to better understand retinal pathology

Heather Day, COA, BS

At the end of the workshop, you'll be able to:

- Describe the basic parts of the posterior ocular anatomy being imaged during OCTs and Fundus Photography.
- List common ocular diseases/conditions for which OCT and Fundus Photography may be ordered.
- Describe 2-4 reasons why an image may be blurred, display artifacts, or be difficult to capture and what to do to maximize image quality for each.
- Demonstrate how to capture clear images on OCT and Fundus Photography.



DID YOU KNOW?

The retina is about the size of a postage stamp.


Inside that small space there are **~3.5 million** retinal pigment epithelium (RPE) cells, **4-5 million** cones, & about **77-107 million** rods!

I make a lot of jokes about eyes...

You could say I have vitreous humor.

What is an OCT?

Optical coherence tomography (OCT) is a non-invasive imaging test that uses light waves to take high-resolution, three-dimensional images of different structures within the eye.



The best way to describe it is like a slice of Birthday Cake!

The OCT allows us to see all of the retinal layers that we wouldn't be able to see without cutting a slice out!

Heidelberg / Zeiss



Heidelberg / Zeiss

What are the differences?


Technical Specs for Heidelberg:

- **Confocal Scanning Laser Ophthalmoscopy:** cSLO provides image detail and clarity not available from fundus photography.
- **TruTrack Active Eye Tracking:** allows for precise, automated retinal follow-up scanning, & retinal thickness measurement reproducibility to 1 micron.
- **AutoRescan:** Automatically places follow-up scans in precisely the same position.
- **SPECTRALIS® Noise Reduction:** Acquires multiple images in the exact same location, differentiating structural information from noise and removes the noise.
- **High resolution OCT images of Retinal layers, Optic Nerve Head, Retina Nerve Fiber Layers, Cornea, Sclera, & Anterior Chamber Angle.**


Technical Specs for Zeiss:

- **FastTrac™:** reduces eye motion artifacts with a proprietary scan acquisition strategy, high speed 20 Hz LSO camera, and single-pass alignment scanning. It also allows you can scan at the highest resolution at the same location at each visit.
- **60-degree OCTA widefield view**
- **High resolution OCT images of Retinal layers, Optic Nerve Head, Retina Nerve Fiber Layers, Cornea, & Anterior Chamber Angle.**

Imaging - Heidelberg




Imaging - Zeiss



What is Fundus Photography?

Fundus photography is a simple, non-invasive imaging test that takes two-dimensional photographs of the back inner wall of your eye, known as the Fundus.



This photo is more like looking at the top of a Birthday Cake.

It allows you to see changes in the anatomy on a more surface level.

Clarus / Optos



Clarus / Optos

What are the differences?

Technical Specs for Clarus:

- Resolution: 7.3 μm
- Field of View:
 - Widefield (one image) - 133°
 - Ultra-widefield (two images) - 200°
 - Montage (up to six images) - up to 267°
- Wavelengths:
 - Red LED: 585 - 640 nm
 - Green LED: 500 - 585 nm
 - Blue LED: 435 - 500 nm
 - Infrared laser diode: 785 nm
- Capture Time: Less than 0.2 seconds


Technical Specs for Optos:

- Resolution: 20 μm
- Field of View:
 - Standard: 200° Single Shot
 - Auto-montage: Up to 220°
- Wavelengths:
 - Red laser: 635 nm
 - Green laser: 532 nm
 - Blue laser: 488 nm
 - Infrared laser: 802 nm
- Capture Time: Less than 0.4 seconds

Image Modalities - Clarus




Ultra-wide True Color




Fundus Autofluorescence




True Color Automontage



Red Channel



Green Channel



Blue Channel


Image Modalities - Optos




Optomap color rgb




Optomap color rg




Optomap Sensory Retina



Optomap choroidal




Optomap green af




Optomap blue af

—
What posterior anatomy do these scans capture?



OCTs / Fundus Photography


- Posterior segment anatomy
 - 10 Retinal Layers
 - Macula
 - Optic Nerve



- Posterior segment anatomy
 - Retina
 - Macula
 - Retinal Blood Vessels
 - Optic Nerve Head
 - Choroid
 - Vitreous




—
What ocular conditions might require an OCT or Fundus Photography?




OCTs / Fundus Photography

Almost all pathology of the posterior segment of the eye can be detected with a fundus camera or an OCT. A few of the most common conditions are:

- macular holes
- Epiretinal Membrane (ERM)
- macular edema
- age-related macular degeneration (AMD)
- geographic atrophy
- central serous retinopathy
- diabetic retinopathy
- vitreous traction
- abnormal blood vessels
- blood vessel blockage
- drusen
- retinal tears
- retinal detachments
- schisis
- lesions and nevi



Why might it be difficult to capture a quality scan?




OCTs / Fundus Photography

OCTs and Fundus Photography rely on light waves to acquire images, therefore they aren't as effective when capturing through conditions that interfere with light passing through the eye; including:


- Poorly dilated pupils
- Clouding or scarring of the cornea
- Dry eyes
- Dense cataracts
- Myopia
- Dense vitreous floaters
- Significant bleeding in the vitreous



Dry Eye



Same eye, 15 minutes later, after PFATs




High Myope

OCTs / Fundus Photography

Additional causes for artifacts or poor images include:

- Poor patient positioning
 - Being too close or too far away
 - Eyelids or Eyelashes in the way
 - The patient's nose casting a shadow
 - Poor pupil alignment
 - Patient movement
- Mask fogging the lens of the instrument
- Dust, smudges, or debris on the camera lens




Mask causing the lens to fog


OCTs / Fundus Photography

How to minimize artifacts


- Ensure proper patient positioning, including sitting straight and aligning the canthus mark with the center of the pupil.
- Instruct the patient to fixate on the target and minimize blinking during the scan, utilizing BSS if necessary.
- Turn the patient's head slightly so their nose is out of the picture.
- Adjust the table height to reduce patient fatigue and movement.
- If needed, hold or tape lids/lashes up and out of the way.
- Properly care for your instrument by keeping it clean and dust free.




Patient Moved




Too Far Away



Lid Obscuring Image



Too Close / Not Centered



Your Turn!

Summing Up



Anatomy

- OCT vs. Fundus Photography
- Retinal Layers
 - Macula
 - Optic Nerve
 - Vitreous



Conditions

- AMD
- ERM
- NPDR
- Drusen
- Macular Holes & Edema



How to Troubleshoot

- Hold/Tape Lids
- Adjust positioning
- Are you too close?
- Are you too far away?
- Are the eyes dry?



How to Capture

- OCTs:
- Heidelberg & Zeiss
- Fundus Photography:
- Optos & Clarus



Any Questions?

Sources

- Bille JF, editor. *High Resolution Imaging in Microscopy and Ophthalmology: New Frontiers in Biomedical Optics* [Internet]. Cham (CH): Springer; 2019. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK554051/> doi:10.1007/978-3-030-16638-0
- Carl Zeiss Meditec, Inc. (2021). ZEISS CLARUS 500. Dublin, California.
- Cleveland Clinic. (2025, July 29). How A picture can help your eyes: What is fundus photography?. Fundus Photography. <https://my.clevelandclinic.org/health/diagnostics/fundus-photography>
- Dong, Zachary & Wollstein, Gadi & Schuman, Joel. (2016). Clinical Utility of Optical Coherence Tomography in Glaucoma. *Investigative Ophthalmology & Visual Science*. 57. OCT556. 10.1167/iovs.16-19933.
- Evranos, Berna & Ugurlu, Nagihan & Faki, Sevgul & Polat, Sefika & Ersoy, Reyhan & Cakir, Bekir. (2021). Retinal layers in prolactinoma patients: a spectral-domain optical coherence tomography study. *International Ophthalmology*. 41. 10.1007/s10792-021-01701-8.
- Heidelberg Engineering Inc. (2025). The DNA of SPECTRALIS. SPECTRALIS OCT - The modular Imaging Platform | Heidelberg Engineering. <https://business-lounge.heidelbergengineering.com/us/en/products/spectralis/spectralis/#product-details>
- Kang, Na & Kim, Hyun Jin & Lee, Jeong-Hee. (2011). The Measurements of Macular Thickness and Volume with SD-OCT in Normal Eyes. *Journal of the Korean Ophthalmological Society*. 52. 1182. 10.3341/jkos.2011.52.10.1182.
- Kemenyova, P & Turcani, Peter & Sutovsky, S & Waczulikova, Iveta. (2014). Optical coherence tomography and its use in optical neuritis and multiple sclerosis. *Bratislavské lekárske listy*. 115. 723-9. 10.4149/BLI_2014_140.
- Liew, A., & Agaian, S. (2025). Comprehensive Survey of OCT-Based Disorders Diagnosis: From Feature Extraction Methods to Robust Security Frameworks. *Bioengineering*, 12(9), 914. <https://doi.org/10.3390/bioengineering12090914>
- Mishra C, Tripathy K. Fundus Camera. [Updated 2023 Aug 25]. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2025 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK558111/>
- Optos. (2023). California. Marlborough, Massachusetts.